BROEN-LAB
FLOW ADJUSTABLE
QUICK CONNECT 5.0

INSTALLATION AND USER MANUAL









BROEN-LAB

Table of contents

Table of Contents	2
Datasheet	3-5
Technical Specifications Built-in Dimensions	3 4-5
Unwrapping the Quick Connect	6-7
070XXX1 (built-in) 070XXX2 (wall/panel mounting, concealed piping) 070XXX3 (wall/panel mounting, exposed piping) 070XXX4 (bench/suspended mounting)	6 6 7 7
Assembly	8-13
Safety Instructions Built-in model Wall bracket Wall/panel mounting, concealed & exposed piping Mounting from frontside without wall bracket Mounting from frontside with wall bracket Bench/suspended mounting	8 9 10 11 11 12 13
Operation	14-15
Trouble Shooting Guide	16
Service	17
Replacement of O-ring on outlet Changing media type	17 17
Spare Part List	18-19

All technical specifications and dimensions are nominal and are subject to change without prior notice.

Datasheet

TECHNICAL SPECIFICATIONS

Standard media¹⁾: Inactive gases

(Air, Vac, O₂, N₂, CO₂, H₂,

Ar & He).

Purity²⁾: 5.0. (99.999%)

Media pressure: -1 to 10 Bar.

-14.5 to 145 psi.

Test pressure: 1.3 x working pressure.

Working

temperature: 0° til 60°C.

32° to 140°F.

Inlet: Female 1/4" NPT.

Visible metal parts: Stainless steel & Chemical


nickel plated brass.

Rubber handle: EPDM.


Metals with

medium contact: Brass.

Gaskets: FKM.



²⁾ Purity 5.0 can only be obtained when the system is purged each time an outlet is connected.




Gas used in test: Compressed Air.

Model used in test: Quick Connect with outlet No. 07 051.00.

 $^{^{\}mbox{\tiny 1)}}$ Other media must first be confirmed by BROEN-LAB .


Datasheet

BUILT-IN Built-in Dimensions



BENCH/SUSPENDED MOUNTING


Built-in Dimensions



Datasheet


WALL/PANEL MOUNTING, CONCEALED PIPING

Built-in Dimensions



WALL/PANEL MOUNTING, EXPOSED PIPING


Built-in Dimensions




BROEN-LAB QUICK CONNECTUnwrapping the Quick Connect

When ordered parts are received, it must be checked that they comply with ordered article. Illustrations below show the contents of delivered articles:

070XXX1 (BUILT-IN)




070XXX2 (WALL/PANEL MOUNTING, CONCEALED PIPING)




BROEN-LAB QUICK CONNECTUnwrapping the Quick Connect

070XXX3 (WALL/PANEL MOUNTING, EXPOSED PIPING)



070XXX4 (BENCH/SUSPENDED MOUNTING)



SAFETY INSTRUCTIONS


Handling of technical gases requires knowledge of the subject, operation instructions and special safety measures. Always check national and local laws, regulations and procedures regarding the installation and use of this type of equipment.

- Changes or modifications are not allowed to be made to the equipment without prior approval from BROFN-I AR
- Improper handling or use can involve risks for the user and other persons as well as damage of the equipment.
- The equipment can only be used for the gas type it
 was manufactured and labelled for. Change can be
 made prior to start of use, provided that the new
 media is listed in the data sheet, or BROEN-LAB has
 approved the new media.
- Observe the maximum pressure & temperature ranges.
- The equipment cannot be used for gases in the liquid phase.


- The equipment must not be exposed to impacts, vibrations or other mechanical forces that can damage the equipment.
- The equipment is designed for use in laboratories.
 The equipment is not designed for installing where there is a possibility of exposure to aggressive media.
- By installation avoid oil and grease contamination (Fire or explosion hazard).
- It is recommended to use a proper leak detection fluid or similar to check for leaks when installing and using the Quick Connect.
- These valves do not contain a relief valve. In case
 of failure, the outlet pressure may rise above the
 adjusted value. Any equipment connected to the
 equipment outlet has to be protected against overpressure.

BUILT-IN MODEL

Built-in dimensions are shown on page 4.




INLET



WALL BRACKET

Wall bracket **07101014** can be used as a tool for marking up the positions of the holes on the panel when models **070XXX2** and **070XXX3** should be mounted from the front or directly on the wall (see pages 11 and 12 for more information).

DESIGN




Assembly

WALL/PANEL MOUNTING, CONCEALED & EXPOSED PIPING

(MOUNTING FROM FRONTSIDE WITHOUT WALL BRACKET)


Built-in dimensions are shown on page 5.



¹⁾ **0700171:** For concealed piping. **0700170:** For exposed piping.

INLET (concealed piping)


INLET (exposed piping)



WALL/PANEL MOUNTING, CONCEALED & EXPOSED PIPING

(MOUNTING FROM FRONTSIDE WITH WALL BRACKET)


Built-in dimensions are shown on page 5.



¹⁾ **0700171:** For concealed piping. **0700170:** For exposed piping.


INLET (concealed piping)

INLET (exposed piping)



BENCH/SUSPENDED MOUNTING

Dimensions are shown on page 4.



Operation

1. Connect outlet and Quick Connect:

- 1.1. Place the outlet in such a way that the first hexagon on the outlet and female hexagon in the Quick Connect fit together. Press the inlet into the Quick Connect until the second hexagon on the outlet touches the Ouick Connect.
- Turn the outlet until the second hexagon on the outlet and the hexagon in the Quick Connect are aligned.
- 1.3. Press the outlet into the Quick Connect.
- 1.4. The outlet is now locked in the Quick Connect.

2. Open valve:

- 2.1. Turn the handle counterclockwise until the desired flow is achieved.
- 2.2. The regulating area is from 0° to 90°.


3. Close valve:


3.1. Turn the handle clockwise until the valve is closed.


4. Release outlet:


- 4.1. Turn the handle clockwise until the outlet is released. The outlet is now parked in the Quick Connect.
- 4.2. Turn the outlet until the first hexagon on the outlet and the hexagon in the Quick Connect are aligned.
- 4.3. Remove the outlet.


Operation




















Trouble shooting guide

ERROR	CAUSE	ACTION
Quick Connect leaks. (Outlet is not plugged in).	1. The Quick Connect leaks between the fitting and inlet.	1 & 2. Retighten the fittings. Apply new sealant to the threads.
	2. The piping leading to the Quick Connect are leaking.	
	3. There is a leak inside the Quick Connect.	3. Replace the Quick Connect with a new one.
The connection between the Quick Connect and the outlet is leaking.	The gasket on the outlet is leaking.	Replace the gasket with a new one.
The outlet is leaking.	There are leaks inside the outlet.	Replace the outlet with a new one.
The handle can be turned without the plug inserted.	The ball lock has been activated by an object other than the outlet.	Turn the handle counter clockwise until the ball lock is released.

Service

REPLACEMENT OF O-RING ON OUTLET

- 1. Remove O-ring as shown.
- Refit a new O-ring. (See spare part list on page 19).



CHANGING THE MEDIA TYPE

In order to be able to change the media type, it is necessary to replace the key and the label.

If the new medium is different from the ones stated on the datasheet, it is necessary to consult BROEN-LAB for confirmation first.

Change of media on a Quick Connect after start of use is not allowed. If the media is changed despite this warning, BROEN-LAB cannot be held responsible for possible accidents!

- 1. Remove label and key. Remove O-ring.
- Refit O-ring on the new key.Refit the new key. Max. Torque: 5Nm.
- 3. Place the new label.
- 4. Replace outlet with new ones suitable for the new media.



Spare part list

Description	Product No. (Similar for all media)
Head cap screw M4x12 (mounting kit)	92E0063
Socket set screw M6x45 (mounting kit)	92E0058
Head cap screw M3x30 (front side mounting)	92E0057
Mounting bracket	0700187
Supporting disk (built-in model)	0210547
Back plate (built-in model)	0700186
Cover (built-in model)	0700173
Cover (exposed piping model)	0700170
Cover (concealed piping model)	0700171

Description	Product No. (Similar for all media)
Mounting column bracket (bench/suspended model)	0210521
Cover for mounting plate (bench/suspended model)	0210522
Mounting plate (bench/suspended model)	0210509
Countersunk screw M4x12 (bench/suspended model)	92G0074
Head cap screw M3x12 (bench/suspended model)	92G0062
Head cap screw M5x20 (bench/suspended model)	92E0064
Head cap screw M5x40 (bench/suspended model)	92E0007
Quick Connect without media key	0700200B
O-ring for media key	92A0183

Description	Product No. (Similar for all media)		
Stainless Steel braided smooth PTEE hose (1.5m) Connections: Male 1/4*NPT - Female 1/4*NPSM Female 1/4*NPSM - Female 1/4*NPSM	93C0735 93C0738		
Swagelok" elbow pipe fitting Material: Stainless Steel	9610113		
Swagelok' Tee pipe fitting Material: Stainless Steel	9610114		
Swagelok' straight pipe fitting Material: Stainless Steel	9610115		
Swagelok' straight tube fitting Material: Stainless Steel Tube diameter: 3 mm 4 mm 6 mm 1/4** 1/8**	96l0100 96l0101 96l0102 96l0103 96l0104		
Swagelok' elbow tube fitting Material: Stainless Steel Tube diameter: 3 mm 4 mm 6 mm 1/4" 1/8"	96l0106 96l0107 96l0108 96l0109 96l0110		

¹⁾ Swagelok – TM Swagelok Company

Spare part list

Description	Label				Media / Pi	roduct No.			
	EN ¹⁾	CA	V	(O ₂)	N ₂	CO	H ₂	Ar	He
	US ²⁾	Air	Vac	O ₂	N_2	CO ₂	H ₂	Ar	He
Label BROEN-LAB	EN ¹⁾ US ²⁾	0700130 0700050	0700131 0700051	0700132 0700052	0700133 0700053	0700134 0700054	0700135 0700055	0700136 0700056	0700137 0700057
Key for quick connect		0700110	0700111	0700112	0700113	0700114	0700115	07000116	0700117
90° outlet, serrated hose nozzle	EN ¹⁾ US ²⁾	0705100 0755100	0705200 0755200	0705300 0755300	0705400 0755400	0705500 0755500	0705600 0755600	0705700 0755700	0705800 0755800
90° outlet, female 1/4″ NPT.	EN ¹⁾ US ²⁾	0705110 0755110	0705210 0755210	0705310 0755310	0705410 0755410	0705510 0755510	0705610 0755610	0705710 0755710	0705810 0755810
90° outlet, female 1/4" NPT.		0707100	0707200	0707300	0707400	0707500	0707600	0707700	0707800
Flowmeter	EN¹) US²)	070513Y ³⁾	-	070533Y ³⁾	070543Y ³⁾ 075543Y ³⁾	070553Y ³⁾	070563Y ³⁾	070573Y ³⁾ 075573Y ³⁾	070583Y ³⁾
Pressure regulator	EN ¹⁾ US ²⁾	070512Z ⁴⁾	-	070532Z ⁴⁾	070542Z ⁴⁾	070552Z ⁴⁾	070562Z ⁴⁾	070572Z ⁴⁾	070582Z ⁴⁾
O-ring for outlet		92A0183	92A0183	92A0182	92A0182	92A0180	92A0180	92A0180	92A0177

¹⁾ Label according to Standard EN13792:2002

²⁾ Label according to American recommendation

³⁾ **Flowmeter:** Y=1: Capacity: 0-5 l/min. Y=2: Capacity: 0-15 l/min. ⁴⁾ **Pressure regulator:** Z=3: 0-500 kPa / 0-3 bar. Z=4: 0-1000 kPa / 0-10bar.

BROEN-LAB A/S

Drejervænget 2 5610 Assens Denmark

Tel: +45 6376 6376 Fax: +45 6471 2476 E-mail: lab@broen-lab.com

www.BROEN-LAB.com

BROEN-LAB